

## PRODUCT DATA SHEET

Standing: 2024-11-15

Stainless steel test tube rack  
with 20 openings Ø18 mm  
Part Number: A000015

### Features

- Stainless steel test tube rack
- suitable for tubes Ø16/17 mm and 15 mL falcons
- Two handles extending from the water surface allow for easy placement and removal.
- Height of the handle: 108 mm
- Distance from bottom: 28 mm



### Technical Features

Overall dimensions (WxDxH)	230 x 47 x 185 mm
Number of samples	20
Material	Stainless steel
Weight	0.4 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser