

PRODUCT DATA SHEET

Standing: 2024-09-17

Stainless steel test tube rack
with 5 openings \varnothing 31 mm
Part Number: A000016

Features

- Stainless steel test tube rack
- Two handles extending from the water surface allow for easy placement and removal.
- Height of the handle: 108 mm
- Distance from bottom: 28 mm



Technical Features

Overall dimensions (WxDxH)	230 x 47 x 185 mm
Number of samples	5
Material	Stainless steel
Weight	0.3 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser