


## PRODUCT DATA SHEET

Standing: 2024-11-04

Stainless steel test tube rack  
with 63 openings Ø31 mm  
Part Number: A000019

### Features

- Stainless steel test tube rack
- Easy positioning and fixing on the four stainless steel pins in the shaker basket.
- Two handles extending from the water surface allow for easy placement and removal.
- Height of the handle: 100 mm
- Distance from bottom: 54 mm



### Technical Features

Overall dimensions (WxDxH)

450 x 300 x 155 mm

Material

Stainless steel

Weight

1.6 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser