

## PRODUCT DATA SHEET

Standing: 2024-11-15

### Shaking tray

Part Number: A000023

#### Features

- Shaking tray made of stainless steel.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and to fix test tube racks A000032 and A000033.
- Perforation: 33 x 21 holes
- Hole grid: 12 mm, hole diameter: 4.2 mm
- Two handles extending from the water surface allow for easy placement and removal.
- When using the tray, the useable depth is reduced by 15 mm: Usable depth: 175 mm



#### Technical Features

Overall dimensions (WxDxH)	450 x 300 x null mm
Usable floor space width	400 mm
Usable floor space depth	256 mm
Material	Stainless steel
Weight	1.0 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser