


## PRODUCT DATA SHEET

Standing: 2024-11-15

Shaking tray for Shakers,  
for all VS 8 models  
Part Number: A000044

### Features

- Shaking tray made of stainless steel.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and separating funnels as well as to fix test tube racks A000059 and A000060.
- Hole grid: 12 mm, hole diameter: 4.2 mm



### Technical Features

Dimensions width	330 mm
Dimensions depth	330 mm
Overall dimensions (WxDxH)	330 x 330 x null mm
Usable floor space width	316 mm
Usable floor space depth	316 mm
Material	Stainless steel

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser