

PRODUCT DATA SHEET

Standing: 2024-09-02

Shaking tray for Shakers, for all VS 15 models and Shaking incubator VS 60 OI
Part Number: A000045

Features

- Shaking tray made of stainless steel.
- With perforated grid to accommodate clamps for Erlenmeyer flasks and separating funnels as well as to fix test tube racks A000059 and A000060.
- Perforation: 33 x 33 holes
- Hole grid: 12 mm, hole diameter: 4.2 mm



Technical Features

Overall dimensions (WxDxH)	450 x 450 x null mm
Usable floor space width	388 mm
Usable floor space depth	388 mm
Material	Stainless steel
Weight	1.3 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Ralf Hermann, Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser