


PRODUCT DATA SHEET

Standing: 2025-01-07

Universal mount for Shakers, for all VS 15 models and Shaking incubator VS 60 OI
Part Number: A000049

Features

- For secure attachment of practically all shaking objects
- Consists of an aluminum attachment with rubber pad, 2 stainless steel longitudinal struts and 6 rubber-coated crossbars
- The longitudinal struts can be mounted in three different height positions (25, 45 and 65 mm above the shelf)
- The mount is placed on the plastic pins of the basic equipment of the shaker



Technical Features

Overall dimensions (WxDxH)	450 x 450 x 80 mm
Usable floor space width	445 mm
Usable floor space depth	400 mm
Weight	4.0 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser