

## PRODUCT DATA SHEET

Standing: 2024-11-15

Test tube rack for max. 24 tubes  
 12 - 17 mm Ø, length 75 - 160 mm  
 Part Number: A000059

### Features

- suitable for test tubes Ø12 - 17 mm, e.g. 15 mL falcon tubes
- The perforated part can be swiveled by 90 ° and is equipped with pressure springs for a tight hold.
- The racks must be screwed onto a tray, e.g. shaking trays A000044 to A000047 for Varioshake. Not suitable for the use in shaking water baths.
- Stainless steel test tube rack



### Technical Features

Overall dimensions (WxDxH)	340 x 60 x 80 mm
Number of samples	24
Weight	0.6 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
 Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
 info@lauda.de • www.lauda.de  
 WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
 Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
 LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
 Sitz Lauda-Königshofen  
 Registergericht Mannheim • HRB 560226

Geschäftsführer:  
 Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
 Dr. Marc Stricker  
 Beirat: Dr. Gerhard Wobser