

PRODUCT DATA SHEET

Standing: 2024-11-15

Test tube rack for max. 16 tubes
 25 - 29 mm Ø, length 75 - 160 mm
 Part Number: A000060

Features

- suitable for test tubes Ø25 - 29 mm, e.g. 50 mL falcon tubes
- The perforated part can be swiveled by 90 ° and is equipped with pressure springs for a tight hold.
- The racks must be screwed onto a tray, e.g. shaking trays A000044 to A000047 for Varioshake. Not suitable for the use in shaking water baths.
- Stainless steel test tube rack



Technical Features

Overall dimensions (WxDxH)	340 x 60 x 80 mm
Number of samples	16
Weight	0.8 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
 Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
 info@lauda.de • www.lauda.de
 WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
 Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
 LAUDA DR. R. WOBSEY Verwaltungs-GmbH
 Sitz Lauda-Königshofen
 Registergericht Mannheim • HRB 560226

Geschäftsführer:
 Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
 Dr. Marc Stricker
 Beirat: Dr. Gerhard Wobser