

## PRODUCT DATA SHEET

Standing: 2024-11-21

Temperature brake for  
data logger LOG100CRYO  
Part Number: A001384

### Features

- Temperature retarder made of aluminum for installation in freezers
- A monitoring sensor can be attached directly in one of the holes in the temperature retarder. This prevents premature alarms when the door/lid is opened briefly
- Three holes  $\varnothing$  3.5mm / 4.2mm / 5mm, each with a hole for a locking screw



### Technical Features

Overall dimensions (WxDxH)

60 x 25 x 25 mm

Weight

0.2 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser