

PRODUCT DATA SHEET

Standing: 2024-09-09

LAUDA Hydro H 6 V

Vaporization water bath 230 V; 50/60 Hz

Part Number: L003091

Features

- Vaporization water bath for gentle steaming work
- Temperature adjustment via rotary knob with temperature scale
- 1-point temperature controller
- Water deficiency protection
- Adjustable water level regulator, sieve bottom and ring set as standard
- Bath interior, cover, sieve plate and heating element are made of stainless steel.
- Multi-piece ring set made of heat-resistant plastic
- Support rods (length 600 mm, diameter 12 mm), made of stainless steel, are included as standard.



Working temperature min.
25 °C



Working temperature max.
100 °C

Technical Features

| | |
|----------------------------|----------------------------------|
| Working temperature range | 25 ... 100 °C |
| Ambient temperature range | 10 ... 40 °C |
| Temperature stability | 3 ± K |
| Heater power max. | 1 kW |
| Power consumption max. | 1 kW |
| Current max. | 16 A |
| Size of bath (Ø) | 131 mm |
| Number of bath openings | 4 |
| Overall dimensions (WxDxH) | 1010 x 270 x 192 mm |
| Weight | 13 kg |
| Power supply | 230 V; 50/60 Hz |
| Power plug | Power cord with plug (NBR 14136) |

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser