

PRODUCT DATA SHEET

Standing: 2025-01-16

Metal tube MTK100 G 3/4" to G3/4"

DN 20, length 1000 mm

Part Number: LZM 075

Features

- Flat sealing

Technical Features

Operating temperature range

-50 ... 150 °C

Connection 1 thread (inside)

G 3/4"

Connection 2 thread (inside)

G 3/4"

Length

100 cm

Material

Stainless steel

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser