

## PRODUCT DATA SHEET

Standing: 2024-11-15

Stainless steel inset rack RN 30  
for 15 tubes, 24-30 mm Dia. 110 mm ID  
Part Number: UE 032

### Features

- Rack with plastic handle(s) for a temperature up to 100 °C
- Distance between levels: 50 mm
- Perforation: 3 x 5 holes



### Technical Features

Overall dimensions (WxDxH)	300 x 145 x 137 mm
Number of samples	15
Immersion depth	110 mm
Material	Stainless steel
Weight	0.8 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser