

PRODUCT DATA SHEET

Standing: 2024-10-31

Stainless steel inset rack RQ 18/2
 for 21 tubes, Ø 14-18 mm 110 mm ID
 Part Number: UG 107

Features

- For comfortable handling in bath thermostats and water baths.
- Rack with plastic handle(s) for a temperature up to 100 °C
- suitable for test tubes with a diameter of 14 - 18 mm
- Distance from bottom: 60 mm



Technical Features

Overall dimensions (WxDxH)	190 x 95 x 203 mm
Number of samples	21
Immersion depth	110 mm
Material	Stainless steel
Weight	0.4 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
 Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
 info@lauda.de • www.lauda.de
 WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
 Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
 LAUDA DR. R. WOBSEY Verwaltungs-GmbH
 Sitz Lauda-Königshofen
 Registergericht Mannheim • HRB 560226

Geschäftsführer:
 Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
 Dr. Marc Stricker
 Beirat: Dr. Gerhard Wobser