

PRODUCT DATA SHEET

Standing: 2024-11-15

Stainless steel inset rack RQ 30
for 10 tubes, Ø 24-30 mm 110 mm ID
Part Number: UG 108

Features

- For comfortable handling in bath thermostats and water baths.
- Rack with plastic handle(s) for a temperature up to 100 °C
- suitable for test tubes with a diameter of 24 - 30 mm
- Distance from bottom: 60 mm

Technical Features

Overall dimensions (WxDxH)	190 x 95 x 203 mm
Number of samples	10
Immersion depth	110 mm
Material	Stainless steel
Weight	0.4 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser